

Sample Prep - Solid Phase Extraction

United Chemical Technologies (UCT) product guide



CLEAN SCREEN®

Drugs of Abuse Columns

- DAU = Acidic, Basic & Neutral Drugs
- THC = Carboxy THC
- THCA = THC Δ9 Carboxylic Acid
- GHB = Gamma-Hydroxybutyrate
- ETG = Ethyl Glucuronide
- BNZ = Benzodiazepine
- CLEAN-THRU® Tips



STYRE SCREEN®

Polymeric Based Columns

- DBX = Copolymeric
- DVB = Divinylbenzene
- BCX = Benzenesulfonic Acid
- C18 = Octadecyl C18
- CCX = Carboxylic Acid
- QAX = Quaternary Amine



CLEAN-UP®

Solid Phase Extraction Columns

- Ion Exchange
- Hydrophobic
- Hydrophilic
- Copolymeric
- Covalent



ENVIRO-CLEAN®

- Products
- Environmental SPE Cartridge
- SPE Inert Glass Syringe Barrels

Sample Prep - Solid Phase Extraction


UCT product guide

First of its kind in the field of forensic and clinical toxicology !

The recently published *Forensic and Clinical Applications of Solid Phase Extraction*, by Michael J. Telepchak, Thomas August and Glynn Chaney, has been met with enthusiasm by those in the SPE field, and has been recommended as a valuable laboratory reference. Dr. Terry Danielson, Ph.D., who reviewed the book for the American Society of Crime Laboratory Directors, calls attention to the "extensive details of many currently available SPE separation procedures", and describes the book as a convenient compendium of SPE technology, and is relevant to the development, implementation and practice of modern SPE appropriate to students, and experienced practitioners."


In the *Canadian Society of Forensic Science Journal*, Dr. Karen Woodall, Ph.D., of the Toronto Centre of Forensic Sciences, calls the book a "must read" for anyone interested in SPE, especially helpful in resolving the day-to-day problems that can occur when using [SPE] and gives some excellent examples of how to deal with some of these occurrences such as recovery variability, contamination, flow, and non-extraction problems."

BOOK "FORENSIC AND CLINICAL APPLICATIONS OF SOLID PHASE EXTRACTION"
P/N : ZZ3801



Reservoirs for Bonded Phase Extractions

Stated Volume (mL)	Tube Configuration	Bed Diameter (mm)	Sorbent Mass (mg)
1	Cylindrical	5.5	50-200
3	Cylindrical	8.5	50-1000
6	Cylindrical	12.5	200-2000
10	Expanded	8.5	50-1000
15	Cylindrical	15.5	500-2000
25	Cylindrical	20	500-5000
75	Cylindrical	27.5	1000-10000
150	Cylindrical	38.0	10000-70000



Sample Prep - Solid Phase Extraction

UCT Clean Screen® columns

Copolymeric bonded phases for drug abuse testing

Analytical demand for more efficient, robust and clean extraction of drugs from biological matrices led to the development of Clean Screen® sorbents. Since 1986, Clean Screen® has led the industry with dependable and reproducible Solid Phase Extraction products and applications. Clean Screen phases are true copolymeric sorbents that contain hydrophobic and ion exchange functional groups uniquely polymerized to a silica substrate. The design and quality of Clean Screen provides superior sample clean up, recovery and reproducibility.

Mixed mode separations allow maximum selectivity for extraction of acids, neutrals and bases. This selectivity makes Clean Screen ideal for both screening and confirmation analysis for virtually all drug categories. Clean Screen DAU, THC, and GHB columns are used extensively by forensic and clinical chemists including :

- Post Mortem Investigations
- Criminal Investigations
- Urine Drug Testing
- Athletic Drug Testing
- Racing Laboratories
- Therapeutic Drug Monitoring
- Medical Drug Screening



Recent additions to this product line include Clean Screen® Ethyl Glucuronide and Clean Screen® Benzodiazepines.

Mechanism of Clean Screen®

Carboxylic acid functionalities present in the sample are ionized. This creates a repulsion between the column and many sample borne interferences, thereby reducing the likelihood of their adsorbing onto the column.

At this pH, ibuprofen & barbiturates are not ionized and are hydrophobically adsorbed onto the column (figure 1). At the same time, drugs with amine functionalities such as cocaine and phencyclidine adsorb onto the column via both hydrophobic and ionic attraction (fig. 1).


The column can then be washed with water or weak aqueous buffers at or below pH 6 without risking loss of the analytes. After drying the column, it is possible to elute the hydrophobically bound analytes using solvents of minimal polarity such as methylene chloride or a hexane/ethyl acetate mixture (fig. 2). Cationic analytes will remain bound to the column. Many compounds of intermediate polarity and potential interferences will also remain on the column. The majority of these potential interferences can be removed by using a methanol wash.




Sample Prep - Solid Phase Extraction

UCT Clean Screen® columns

Cationic analytes bound to the column can be eluted after another drying step. The drying steps are necessary to remove water which would have prevented the water-immiscible elution solvents from optimally interacting with the analytes (fig. 3).



To elute the cationic analytes, an organic solution with a high pH (between 11 & 12) should be used. A methylene chloride isopropanol-ammonium hydroxide mixture will simultaneously disrupt these ionic interactions and successfully elute the desired compound (fig. 4).



Clean Screen® DAU

This column is copolymerized on a rigid, purified silica gel support. The two functional groups include a reverse phase, and an ion exchanger, benzenesulfonic acid. This column is commonly used for analyzing a wide range of drugs of abuse, including acidic, basic & neutral drugs.

Application : Dual functionality for weak bases and hydrophobic compounds.

P/N	Weight /Vol.	Qty
CSDAU131	130 mg/1 mL	100
CSDAU133	130 mg/3 mL	50
CSDAU203	200 mg/3 mL	50
CSDAU303	300 mg/3 mL	50
CSDAU503	500 mg/3 mL	50
CSDAU206	200 mg/6 mL	50
CSDAU506	500 mg/6 mL	50
CSDAU1M6	1 g/6 mL	30
ZSDAU005	50 mg/10 mL	50
ZSDAU013	130 mg/10 mL	50
ZSDAU020	200 mg/10 mL	50
CSDAU515	500 mg/15 mL	50

Clean Screen® THC

This column is copolymerized on a rigid, purified silica gel support. The two functional groups include a reverse phase, and an ion exchanger, primary amine. This column is commonly used for analyzing THC and its metabolites.


Application : Dual functionality for acids and hydrophobic compounds.

P/N	Weight /Vol.	Qty
CSTHC131	130 mg/1 mL	100
CSTHC203	200 mg/3 mL	50
CSTHC303	300 mg/3 mL	50
CSTHC503	500 mg/3 mL	50
CSTHC206	200 mg/6 mL	50
CSTHC506	500 mg/6 mL	50
CSTHC1M6	1 g/6 mL	30
ZSTHC013	130 mg/10 mL	50
ZSTHC020	200 mg/10 mL	50
CSTHC515	500 mg/15 mL	50

Clean Screen® GHB

The small polar nature of the molecule and the lack of a UV chromophore complicate the chromatographic and spectrophotometric analysis of GHB. Chemically, GHB is unstable and readily forms Gamma-butyrolactone when heated in acid conditions. Most analytical methods are based upon the interconversion to the lactone and chemical derivatization to form the TMS derivative. This column is for the extraction of free GHB.

P/N	Weight /Vol.	Qty
CSGHB203	200 mg/3 mL	50
ZSGHB020	200 mg/10 mL	50
ZCGHB020	200 mg/10 mL	50



Sample Prep - Solid Phase Extraction

UCT Clean-Up® columns

Hydrophobic Extraction Columns



This sorbent is composed of a silica backbone bonded with hydrocarbon chains. It is used to extract compounds which exhibit non-polar or neutral characteristics out of complex matrices. The C18 phase is the most widely used for non-polar interactions because of its nonselective nature; C18 will extract a large number of compounds with differing chemical properties. To enhance selectivity, UCT offers a wide range of hydrophobic sorbents, from C2 to C20. Multiple chain configurations are available for some sorbents. Endcapped or unendcapped sorbents are available for all chain lengths.

Mechanism of hydrophobic bonding


Compounds are retained by non-polar interactions from polar solvents or matrix environments. They are bound by dispersion forces / van der Waals forces. Elution, or disruption of the non-polar interactions is achieved by solvents or solvent mixtures with sufficient non-polar character. Some polar solvents, such as acetonitrile have enough non-polar characteristics to disrupt non-polar binding to cause elution of a compound from the sorbent. Methanol can be used as well, although it should be noted that it will take off both polar & non-polar analytes of interest & interferences.

Unendcapped vs. Endcapped

Bonded phases are manufactured by the reaction of organosilanes with activated silica. During the polymerization reaction of carbon chains to the silica backbone, a very stable silyl ether linkage forms. Our unendcapped columns allow hydroxyl sites to remain, thus making these columns slightly hydrophilic. In order to decrease this slight polarity, these hydroxyl sites are deactivated. Proprietary bonding techniques ensure that these sites are 100% reacted, leading to a complete endcapping. Because there are no hydroxyl sites left, our endcapped columns are more hydrophobic than our unendcapped columns.



Example of Hydrophobic Bonding



Functionalized hydrophobic silica based phases

Sorbent	Product code	Structure	% Carbon
C2 ethyl	C02	-SiCH ₂ CH ₃	6.60
C3 propyl	C03	-Si-(CH ₂) ₂ CH ₃	7.60
C4 n-butyl	Cn4	-Si-(CH ₂) ₃ CH ₃	8.50
Ci4 isobutyl	Ci4	-Si-CH ₂ CH(CH ₃) ₂	8.80
Ct4 tertiary butyl	Ct4	-Si-C(CH ₃) ₃	8.50
C5 pentyl	C05	-Si-(CH ₂) ₄ CH ₃	9.50
C6 hexyl	C06	-Si-(CH ₂) ₅ CH ₃	11.00
C7 heptyl	C07	-Si-(CH ₂) ₆ CH ₃	11.00
C8 octyl	C08	-Si-(CH ₂) ₇ CH ₃	11.10
C10 decyl	C10	-Si-(CH ₂) ₉ CH ₃	15.70
C12 dodecyl	C12	-Si-(CH ₂) ₁₁ CH ₃	not tested
C18 octadecyl	C18	-Si-(CH ₂) ₁₇ CH ₃	21.70
C20 eicosyl	C20	-Si-(CH ₂) ₁₉ CH ₃	24.30
C30 tricontyl	C30	-Si-(CH ₂) ₂₉ CH ₃	26.00
Cyclohexyl	CYH1	-Si-(CH)	11.60
Phenyl	PHY1	-Si-(PH)	11.00

Sample Prep - Solid Phase Extraction

UCT Clean-Up® columns

Weight/Vol.	Sorbent	Endcapped	Unendcapped	Sorbent	Endcapped	Unendcapped	Sorbent	Endcapped	Unendcapped	Qty
50mg/1mL	C2, Ethyl	CEC021L1	CUC021L1	C4, n-Butyl*	CECN41L1	CUCN41L1	C6, Hexyl	CEC061L1	CUC061L1	100
100mg/1mL		CEC02111	CUC02111		CECN4111	CUCN4111		CEC06111	CUC06111	100
100mg/3mL		CEC02113	CUC02113		CECN4113	CUCN4113		CEC06113	CUC06113	50
200mg/3mL		CEC02123	CUC02123		CECN4123	CUCN4123		CEC06123	CUC06123	50
500mg/3mL		CEC02153	CUC02153		CECN4153	CUCN4153		CEC06153	CUC06153	50
500mg/6mL		CEC02156	CUC02156		CECN4156	CUCN4156		CEC06156	CUC06156	50
1g/6mL		CEC021M6	CUC021M6		CECN41M6	CUCN41M6		CEC061M6	CUC061M6	30
100mg/10mL		CEC0211Z	CUC0211Z		CECN411Z	CUCN411Z		CEC0611Z	CUC0611Z	50
200mg/10mL		CEC0212Z	CUC0212Z		CECN412Z	CUCN412Z		CEC0612Z	CUC0612Z	50
500mg/10mL		CEC0215Z	CUC0215Z		CECN415Z	CUCN415Z		CEC0615Z	CUC0615Z	50
50mg/1mL	C3, Propyl	CECN31L1	CUCN31L1	C5, Pentyl	CEC051L1	CUC051L1	C7, Heptyl	CEC071L1	CUC071L1	100
100mg/1mL		CECN3111	CUCN3111		CEC05111	CUC05111		CEC07111	CUC07111	100
100mg/3mL		CECN3113	CUCN3113		CEC05113	CUC05113		CEC07113	CUC07113	50
200mg/3mL		CECN3123	CUCN3123		CEC05123	CUC05123		CEC07123	CUC07123	50
500mg/3mL		CECN3153	CUCN3153		CEC05153	CUC05153		CEC07153	CUC07153	50
500mg/6mL		CECN3156	CUCN3156		CEC05156	CUC05156		CEC07156	CUC07156	50
1g/6mL		CECN31M6	CUCN31M6		CEC051M6	CUC051M6		CEC071M6	CUC071M6	30
100mg/10mL		CECN311Z	CUCN311Z		CEC0511Z	CUC0511Z		CEC0711Z	CUC0711Z	50
200mg/10mL		CECN312Z	CUCN312Z		CEC0512Z	CUC0512Z		CEC0712Z	CUC0712Z	50
500mg/10mL		CECN315Z	CUCN315Z		CEC0515Z	CUC0515Z		CEC0715Z	CUC0715Z	50

*Available on request C4 Isobutyl and C4 Tertiary Butyl

Weight/Vol.	Sorbent	Endcapped	Unendcapped	Sorbent	Endcapped	Unendcapped	Sorbent	Endcapped	Unendcapped	Qty
50mg/1mL	C8, Octyl	CEC081L1	CUC081L1	C12, nDodecyl	CEC121L1	CUC121L1	Cyclohexyl	CECYH1L1	CUCYH1L1	100
100mg/1mL		CEC08111	CUC08111		CEC12111	CUC12111		CECYH111	CUCYH111	100
100mg/3mL		CEC08113	CUC08113		CEC12113	CUC12113		CECYH113	CUCYH113	50
200mg/3mL		CEC08123	CUC08123		CEC12123	CUC12123		CECYH123	CUCYH123	50
500mg/3mL		CEC08153	CUC08153		CEC12153	CUC12153		CECYH153	CUCYH153	50
500mg/6mL		CEC08156	CUC08156		CEC12156	CUC12156		CECYH156	CUCYH156	50
1g/6mL		CEC081M6	CUC081M6		CEC121M6	CUC121M6		CECYH1M6	CUCYH1M6	30
100mg/10mL		CEC0811Z	CUC0811Z		CEC1211Z	CUC1211Z		CECYH11Z	CUCYH11Z	50
200mg/10mL		CEC0812Z	CUC0812Z		CEC1212Z	CUC1212Z		CECYH12Z	CUCYH12Z	50
500mg/10mL		CEC0815Z	CUC0815Z		CEC1215Z	CUC1215Z		CECYH15Z	CUCYH15Z	50
50mg/1mL	C10, nDecyl	CEC101L1	CUC101L1	C18, Octadecyl	CEC181L1	CUC181L1	Phenyl	CEPHY1L1	CUPHY1L1	100
100mg/1mL		CEC10111	CUC10111		CEC18111	CUC18111		CEPHY111	CUPHY111	100
100mg/3mL		CEC10113	CUC10113		CEC18113	CUC18113		n.a.	n.a.	50
200mg/3mL		CEC10123	CUC10123		CEC18123	CUC18123		CEPHY123	CUPHY123	50
500mg/3mL		CEC10153	CUC10153		CEC18153	CUC18153		CEPHY153	CUPHY153	50
500mg/6mL		CEC10156	CUC10156		CEC18156	CUC18156		CEPHY156	CUPHY156	50
1g/6mL		CEC101M6	CUC101M6		CEC181M6	CUC181M6		CEPHY1M6	CUPHY1M6	30
100mg/10mL		CEC1011Z	CUC1011Z		CEC1811Z	CUC1811Z		CEPHY11Z	CUPHY11Z	50
200mg/10mL		CEC1012Z	CUC1012Z		CEC1812Z	CUC1812Z		CEPHY12Z	CUPHY12Z	50
500mg/10mL		CEC1015Z	CUC1015Z		CEC1815Z	CUC1815Z		CEPHY15Z	CUPHY15Z	50

Sample Prep - Solid Phase Extraction

UCT Clean-Up® columns

Hydrophilic Normal Phase Columns

This sorbent is composed of a silica backbone bonded with carbon chains containing polar functional groups. Groups which will possess such polarity include amines, hydroxyls and carbonyls.


Functionalized hydrophilic silica based phases

Sorbent	Product code	Structure	% Carbon
Silica	SIL1	-SiOH	N/A
Diol	DOL1	-Si-(CH ₂) ₃ OCH ₂ CHOHCH ₂ OH	8.00
Cyanopropyl	CNP1	-Si-(CH ₂) ₃ CN	6.90
Florisil®	FLS		N/A
Alumina, Acidic	ALA		N/A
Alumina, Neutral	ALN		N/A
Alumina, Basic	ALB		N/A
Carbon	CARB		N/A

Mechanism of hydrophilic bonding

Compounds are retained on hydrophilic sorbents through polar interactions including hydrogen bonding, pi-pi or dipole-dipole interaction. These types of interactions occur when a distribution of electrons between individual atoms in functional groups is unequal, causing negative and positive polarity. Compounds typically extracted on a hydrophilic column include analytes which have polar groups, including amines, hydroxyls and carbonyls. Elution is performed by strong polar solvents.

Example of Hydrophilic Bonding




Weight /Vol.	Sorbent	P/N	Sorbent	P/N	Sorbent	P/N	Sorbent	P/N	Qty
50mg/1mL	Unbonded Silica	CUSIL1L1	Florisil®	CUFLS1L1	Alumina, Acidic	CUALA1L1	Alumina, Basic	CUALB1L1	100
100mg/1mL		CUSIL111		CUFLS111		CUALA111		CUALB111	100
200mg/3mL		CUSIL123		CUFLS123		CUALA123		CUALB123	50
500mg/3mL		CUSIL153		CUFLS153		CUALA153		CUALB153	50
500mg/6mL		CUSIL156		CUFLS156		CUALA156		CUALB156	50
1g/6mL		CUSIL1M6		CUFLS1M6		CUALA1M6		CUALB1M6	30
100mg/10mL		CUSIL11Z		CUFLS11Z		CUALA11Z		CUALB11Z	50
200mg/10mL		CUSIL12Z		CUFLS12Z		CUALA12Z		CUALB12Z	50
500mg/10mL		CUSIL15Z		CUFLS15Z		CUALA15Z		CUALB15Z	50
50mg/1mL	Alumina, neutral	CUALN1L1	CN, Cyanopropyl Endcapped	CECNP1L1	CN, Cyanopropyl Unendcapped	CUCNP1L1	Diol	CUDOL1L1	100
100mg/1mL		CUALN111		CECNP111		CUCNP111		CUDOL111	100
200mg/3mL		CUALN123		CECNP123		CUCNP123		CUDOL123	50
500mg/3mL		CUALN153		CECNP153		CUCNP153		CUDOL153	50
500mg/6mL		CUALN156		CECNP156		CUCNP156		CUDOL156	50
1g/6mL		CUALN1M6		CECNP1M6		CUCNP1M6		CUDOL1M6	30
100mg/10mL		CUALN11Z		CECNP11Z		CUCNP11Z		CUDOL11Z	50
200mg/10mL		CUALN12Z		CECNP12Z		CUCNP12Z		CUDOL12Z	50
500mg/10mL		CUALN15Z		CECNP15Z		CUCNP15Z		CUDOL15Z	50

Sample Prep - Solid Phase Extraction

UCT Clean-Up® columns

Ion Exchange extraction columns

This sorbent is composed of a silica backbone bonded with a carbon chain terminated by a negatively or positively charged functional group. Ion exchange interactions occur between a sorbent that carries a charge and a compound of opposite charge. This electrostatic interaction is reversible by neutralizing the sorbent and/or analyte. Ion exchange bonds can also be disrupted by introduction of a "counter ion" to compete with the analyte for binding sites on the sorbent.



Sorbent	Code	Structure	Pka	% Carbon	meq /g
Aminopropyl (1° amine)	NAX1	-Si-(CH ₂) ₃ NH ₂	9.8	6.65	0.310
N-2 Aminoethyl (1° & 2° amine)	PSA1	-Si-(CH ₂) ₃ NH(CH ₂) ₂ NH ₂	10.1, 10.9	9.70	0.320
Diethylamino (3° amine)	DAX1	-Si-(CH ₂) ₃ N(CH ₂ CH ₃) ₂	10.6	8.40	0.280
Quaternary Amine Chloride	QAX1	-Si-(CH ₂) ₃ N+(CH ₃) ₃ Cl ⁻	always charged	8.40	0.250
Quaternary Amine Hydroxide	CHQAX1	-Si-(CH ₂) ₃ N+(CH ₃) ₃ CH ₃ CO ₂ ⁻	always charged	8.40	0.250
Quaternary Amine Acetate	CAQAX1	-Si-(CH ₂) ₃ N+(CH ₃) ₃ OH ⁻	always charged	8.40	0.250
Quaternary Amine Formate	CFQAX1	-Si-(CH ₂) ₃ N+(CH ₃) ₃ CHO ₂ ⁻	always charged	8.40	0.250
Cation Exchange					
Carboxylic Acid	CCX1	-Si-CH ₂ COOH	4.8	9.10	0.170
Propylsulfonic Acid	PCX1	-Si-(CH ₂) ₃ SO ₃ H	<1	7.10	0.180
Benzenesulfonic Acid	BCX1	-Si-(CH ₂) ₂ -(PH)-SO ₃ H	always charged	11.00	0.320
Benzenesulfonic Acid High Load	BCXH1	-Si-(CH ₂) ₂ -(PH)-SO ₃ H	always charged	15.00	0.650

Mechanism of Ion Exchange bonding

Compounds are retained on the sorbent through ionic bonds. Therefore, it is essential that the sorbent and the analyte to be extracted are charged. Generally, the number of molecules with charged cationic groups increases at pH values below the molecule's pKa value. The number of molecules with charged anionic groups decreases at pH values below the molecule's pKa value. To ensure 99% or more ionization, the pH should be at least two pH units below the pKa of the cation and two pH units above the pKa of the anion. Elution occurs by using a solvent to raise the pH above the pKa of the cationic group or to lower the pH below the pKa of the anion to disrupt retention. At this point, the sorbent or compound will be neutralized.

Weight /Vol.	Sorbent	P/N	Sorbent	P/N	Sorbent	P/N	Sorbent	P/N	Qty
50mg/1mL	Aminopropyl	CUNAX1L1	"PSA N-2 Aminoethyl"	CUPSA1L1	Diethylamino	CUDAX1L1	"Quaternary Amine with Chloride Counter Ion**"	CUQAX1L1	100
100mg/1mL		CUNAX111		CUPSA111		CUDAX111		CUQAX111	100
200mg/3mL		CUNAX123		CUPSA123		CUDAX123		CUQAX123	50
500mg/3mL		CUNAX153		CUPSA153		CUDAX153		CUQAX153	50
500mg/6mL		CUNAX156		CUPSA156		CUDAX156		CUQAX156	50
1g/6mL		CUNAX1M6		CUPSA1M6		CUDAX1M6		CUQAX1M6	30
100mg/10mL		CUNAX11Z		CUPSA11Z		CUDAX11Z		CUQAX11Z	50
200mg/10mL	Carboxylic Acid	CUNAX12Z	Propyl sulfonic Acid	CUPSA12Z	Benzene- -sulfonic Acid	CUDAX12Z	Benzene- -sulfonic Acid High Load	CUQAX12Z	50
500mg/10mL		CUNAX15Z		CUPSA15Z		CUDAX15Z		CUQAX15Z	50
50mg/1mL		CUCCX1L1		CUPCX1L1		CUBCX1L1		CUBCX1H1L1	100
100mg/1mL		CUCCX111		CUPCX111		CUBCX111		CUBCX1HL11	100
200mg/3mL		CUCCX123		CUPCX123		CUBCX123		CUBCX1HL23	50
500mg/3mL		CUCCX153		CUPCX153		CUBCX153		CUBCX1HL53	50
500mg/6mL		CUCCX156		CUPCX156		CUBCX156		CUBCX1HL56	50
1g/6mL	Carboxylic Acid	CUCCX1M6	Propyl sulfonic Acid	CUPCX1M6	Benzene- -sulfonic Acid	CUBCX1M6	Benzene- -sulfonic Acid High Load	CUBCX1HLM6	30
100mg/10mL		CUCCX11Z		CUPCX11Z		CUBCX11Z		CUBCX1HL1Z	50
200mg/10mL		CUCCX12Z		CUPCX12Z		CUBCX12Z		CUBCX1HL2Z	50
500mg/10mL		CUCCX15Z		CUPCX15Z		CUBCX15Z		CUBCX1HL5Z	50

**Available with Acetate , Hydroxide & Formate Counter Ion

Sample Prep - Solid Phase Extraction

UCT Clean-Up® columns

Copolymeric Extraction columns (ion Exchange + hydrophobicity)


This sorbent is composed of a silica backbone with two types of functional chains attached - an ion exchanger or polar chain and a hydrophobic carbon chain. Our copolymeric phases are produced in a way to allow for equal parts of each functional group to attach to the silica backbone. This copolymerization technique yields reproducible bonded phases and unique copolymeric chemistries which allow the controlled use of mixed mode separation mechanisms. This type of dual chemistry is beneficial especially when one is looking for both a neutral & charged compound. This is common when a neutral parent drug metabolizes and becomes a charged compound.



Functionalized mixed mode silica based phases



Sorbent	Product code	Structure	% Carbon	meq/g
Aminopropyl + C8	NAX2	-Si-(CH ₂) ₃ NH ₂ & -Si-(CH ₂) ₇ CH ₃	12.3	0.163
Quaternary Amine + C8	QAX2	-Si-(CH ₂) ₃ N ⁺ (CH ₃) ₃ & -Si-(CH ₂) ₇ CH ₃	13.60	0.160
Carboxylic Acid + C8	CCX2	-Si-CH ₂ COOH & -Si-(CH ₂) ₇ CH ₃	12.50	0.105
Propylsulfonic Acid + C8	PCX2	-Si-(CH ₂) ₃ SO ₃ H & -Si-(CH ₂) ₇ CH ₃	14.62	0.114
Benzenesulfonic Acid + C8	BCX2	-Si-(CH ₂) ₂ -(PH)-SO ₃ H & -Si-(CH ₂) ₇ CH ₃	12.30	0.072
Cyanopropyl + C8	CNP2	-Si-(CH ₂) ₃ CN & -Si-(CH ₂) ₇ CH ₃	14.60	0.163
Cyclohexyl + C8	CYH2	-Si-(PH) & -Si-(CH ₂) ₇ CH ₃	N/A	N/A

Example of Copolymeric Bonding



Mechanism of mixed mode bonding

Using a sample composed of a theoretical neutral parent drug and its charged (acidic) metabolite, it is applied at a pH of 6 (figure 1). At this pH, many amine groups are positively charged. Since the column is also positively charged, compounds with this chemistry (cations) are repelled. Depending on the pKa of the metabolite, carboxylic acid groups may be negatively charged, allowing the metabolite to bond to the positively charged sorbent. Since the column also possesses a hydrophobic chain, the neutral parent drug also bonds to the column. Water or a weak aqueous buffer (pH6) washes away hydrophilically bound interferences (figure 2). The column is then dried, careful to free the column of any residual aqueous phase that would interfere with elution.




Sample Prep - Solid Phase Extraction

UCT Clean-Up® columns


Elution 1
(Figure 3)

The hydrophobically bound neutral parent drug is eluted with a solvent of minimal polarity, such as hexane/ethyl acetate - 80:20.



Elution 2
(Figure 4)

The final elution employs an acid to neutralize the charge of acidic analytes. Ionic interaction is released, and analytes are eluted in an appropriate solvent mixture.



Weight /Vol.	Sorbent	P/N	Sorbent	P/N	Sorbent	P/N	Sorbent	P/N	Qty
50 mg/1 mL	Hydrophobic	CUPCX2L1	Hydrophobic	CUCCX2L1	Hydrophobic	CUBCX2L1	Octadecyl	CUBCX3L1	100
100 mg/1 mL	plus	CUPCX211	plus	CUCCX211	plus	CUBCX211	plus	CUBCX311	100
200 mg/3 mL	Propylsulfonic	CUPCX223	Carboxylic	CUCCX223	Benzene-	CUBCX223	Benzene-	CUBCX323	50
500 mg/3 mL	Acid	CUPCX253	Acid	CUCCX253	-sulfonic	CUBCX253	-sulfonic	CUBCX353	50
500 mg/6 mL		CUPCX256		CUCCX256	Acid	CUBCX256	Acid	CUBCX356	50
1 g/6 mL		CUPCX2M6		CUCCX2M6		CUBCX2M6		CUBCX3M6	30
100 mg/10 mL		CUPCX21Z		CUCCX21Z		CUBCX21Z		CUBCX31Z	50
200 mg/10 mL		CUPCX22Z		CUCCX22Z		CUBCX22Z		CUBCX32Z	50
500 mg/10 mL		CUPCX25Z		CUCCX25Z		CUBCX25Z		CUBCX35Z	50
50 mg/1 mL	Hydrophobic	CUPSA2L1	Octadecyl	CUPSA3L1	Hydrophobic	CUQAX2L1	Hydrophobic	CUNAX2L1	100
100 mg/1 mL	plus	CUPSA211	plus	CUPSA311	plus	CUQAX211	plus	CUNAX211	100
200 mg/3 mL	N-2	CUPSA223	N-2	CUPSA323	Quaternary	CUQAX223	Aminopropyl	CUNAX223	50
500 mg/3 mL	Aminoethyl	CUPSA253	Aminoethyl	CUPSA353	Amine	CUQAX253		CUNAX253	50
500 mg/6 mL		CUPSA256		CUPSA356		CUQAX256		CUNAX256	50
1 g/6 mL		CUPSA2M6		CUPSA3M6		CUQAX2M6		CUNAX2M6	30
100 mg/10 mL		CUPSA21Z		CUPSA31Z		CUQAX21Z		CUNAX21Z	50
200 mg/10 mL		CUPSA22Z		CUPSA32Z		CUQAX22Z		CUNAX22Z	50
500 mg/10 mL		CUPSA25Z		CUPSA35Z		CUQAX25Z		CUNAX25Z	50

Sample Prep - Solid Phase Extraction

UCT Styre Screen® Polymeric Resin columns

Styre Screen® extraction columns contain an ultra clean, highly cross-linked styrene and divinylbenzene copolymer sorbent that is functionalized with both a reverse phase, hydrophobic component and a strong cation exchanger. High & reproducible recoveries for acidic, neutral and basic compounds are achievable with a single column. The Styre Screen® particles have an average particle size of 30 microns and a very high analyte capacity making them ideal for standard solid phase extraction applications. The increased analyte capacity means that less sorbent bed mass is needed which results in faster flow rates and less solvent use. Higher throughput and less solvent waste disposal translate into significant savings in both time and money. In addition, no conditioning steps are required for most drugs of abuse applications.



Advantage

- No conditioning steps
- Copolymer allows for extraction of acids, neutrals and bases
- High and reproducible recoveries
- Clean extractions
- Highly cross-linked styrene/divinylbenzene polymer
- Reduction in sorbent mass
- Faster flow rates
- pH stable (1 to 14)
- Reduction in solvent use
- High sorbent capacity
- Methods for NIDA/SAMHSA 5 Drugs

Weight /Vol.	Sorbent	P/N	Application	Qty
30 mg/1 mL	DBX - Benzenesulfonic Acid + C18	SSDBX031	Dual functionality for weak acids and hydrophobic compounds.	100
30 mg/3 mL		SSDBX033		50
50 mg/6 mL		SSDBX056		50
30 mg/1 mL	DVB - Polystyrene Divinylbenzene	SSDBX031	n.a.	100
30 mg/3 mL		SSDBX033	n.a.	50
50 mg/6 mL		SSDBX056	n.a.	50
30 mg/1 mL	BCX - Reverse Phase	SSBCX031	Scavenger for amines, alcohols and other nucleophiles.	100
30 mg/3 mL		SSBCX033		50
50 mg/6 mL		SSBCX056		50
30 mg/1 mL	C18 -Reverse Phase	SSC18031	Removes hydrophobic impurities, de-salting and purification of hydrophobic compounds.	100
30 mg/3 mL		SSC18033		50
50 mg/6 mL		SSC18056		50
30 mg/1 mL	CCX - Carboxylic Acid	SSCCX031	Scavenger for strong anions (Quaternary amines or metals)	100
30 mg/3 mL		SSCCX033		50
50 mg/6 mL		SSCCX056		50
30 mg/1 mL	QAX - Quaternary Amine	SSQAX031	"Removes large or more hydrophobic compounds."	100
30 mg/3 mL		SSQAX033		50
50 mg/6 mL		SSQAX056		50



Sample Prep - Solid Phase Extraction

UCT Enviro-Clean® - Universal cartridges



The Enviro-Clean® Universal Cartridge is the choice of modern contract labs. This inexpensive, easy to use cartridge provides consistent extractions with clean blanks. Built in flow control allows for consistent flow rates. Enviro-Clean® sorbents UCT polypropylene, and PTFE frits offer a clean blank with every batch. Designed for the environmental lab, the cartridge is made to handle large volumes of waste water. An optional bottle holder is available for continuous feed from Boston Round and wide mouth bottles.



SPE-DEX® is a registered trademark of Horizon Technology, Inc.

Product Name	P/N	Weight /Vol.	Description / Application	Qty
UNIVERSAL C18	ECUNIC18	1100 mg/83 mL	1100 mg of endcapped C18 for pesticides, PCBs and a large assortment of applications.	8
UNIVERSAL 525	ECUNI525	1500 mg/83 mL	1500 mg of our special C18 blend. This cartridge is specifically designed for EPA Method 525.	8
UNIVERSAL PAH / DRO	ECUNIPAH	2000 mg/83 mL	2000 mg of C18 specifically designed for PAH extraction.	8
UNIVERSAL OIL & GREASE	ECUNIOAG	4000 mg/83 mL	4000 mg of large particle C18 with an assortment of PE frit filters. No more liquid/liquid emulsions or clogged disks.	15

The cartridge will fit all standard manifolds and disk manifolds with adapter.

Sample Prep - Solid Phase Extraction

UCT dispersive SPE – “QuEChERS”

Quick, Easy, Cheap, Effective, Rugged and Safe Approach for determining pesticide residues in fruits, vegetables and other foods.

The QuEChERS method is gaining in popularity around the world as the method of choice for food testing. The QuEChERS method offers the advantages of high recoveries, accurate results, high sample throughput, low solvent and glassware usage (no chlorinated solvents), less labor and bench space and lower reagent costs. Organic acids and other potential contaminants are removed during the cleanup process.

UCT provides a variety of QuEChERS products containing primary secondary amine (PSA), C18, magnesium sulfate anhydrous and graphitized carbon black. These products are used in the method's clean-up step. Bulk, pre-cleaned magnesium sulfate anhydrous is available for the extraction part of the method.

PSA is used to remove various polar organic acids, polar pigments, some sugars and fatty acid co-extractables from QuEChERS extracts. Combined with C18, samples containing less than 1% lipids can be cleaned of most lipids and sterols. Graphitized carbon is used to remove sterols and pigments such as chlorophyll. The downside to carbon is its ability to retain planar molecules. Schenck and Vega (April 2001) reported that 3/1 acetone : toluene performed well at eluting many compounds from carbon.



P/N	Description	Qty
CUMPSCB2CT	2 mL micro-centrifuge tubes with 150 mg Anhydrous Magnesium Sulfate, 50 mg PSA & 50mg Carbon	100
CUMPS2CT	2 mL micro-centrifuge tubes with 150 mg Anhydrous Magnesium Sulfate, 50 mg PSA	100
CUMPS18CT	2 mL micro-centrifuge tubes with 150 mg Anhydrous Magnesium Sulfate, 50 mg PSA & 50 mg endcapped C18	100
ECMPSCB15CT	15 mL centrifuge tubes with 900 mg Anhydrous Magnesium Sulfate, 300 mg PSA & 50 mg endcapped C18	50
ECMPSC1815CT	15 mL centrifuge tubes with 900 mg Anhydrous Magnesium Sulfate, 300 mg PSA & 150 mg endcapped C18	50
ECPSACB6	6 mL columns with 400 mg PSA on bottom, 200 mg Graphitized Carbon-Black on top, separated by a Teflon frit*	30
ECPSACB256	6 mL columns with 250 mg Graphitized Carbon on top, 500 mg PSA on the bottom, separated with a Teflon frit*	30
ECPSACB506	6 mL columns with 500 mg Graphitized Carbon on top, 500 mg PSA on the bottom, separated with a Teflon frit*	30
ECMSSA50CT	50 mL PP centrifuge tube with 6 g Anhydrous Magnesium Sulfate, 1.5 g Anhydrous Sodium Acetate	250
ECMSSC50CT	50 mL PP centrifuge tube with 4 g Anhydrous Magnesium Sulfate, 1 g NaCl	250

Example of procedure

1. Transfer 15 g of homogenized sample into a 50 ml FEP centrifuge tube.
2. Add 15 ml of 1% acetic acid in acetonitrile, 1.5 g sodium acetate anhydrous, 6 g of UCT magnesium sulfate anhydrous and an internal standard.
3. Shake vigorously for 1 minute.
4. Centrifuge for 3 minutes at 3700 rpm.
5. Transfer an aliquot of the supernatant to the UCT product.
6. Shake for 1 minute.
7. Centrifuge for 3 minutes at 3700 rpm.
8. Analyze.

* Products available with Polyethylene or Teflon frits. Choice depends application and price requirements.

Sample Prep - Solid Phase Extraction

UCT Enviro-Clean® columns

Polypropylene and Inert Glass Extraction columns

Enviro-Clean® solid phase extraction columns are designed especially for the isolation and separation of environmental analytes such as pesticides, herbicides, polyaromatic hydrocarbons, polychlorinated biphenyls and other environmentally related compounds.

Enviro-Clean® offers a selection of high quality solid phase extraction columns geared to support the environmental chemist with a very broad range of analytical applications. The most important function of the solid phase extraction column for the environmental chemist is the clean separation of an analyte from a variety of compounds. An important function of the extraction column is that it will concentrate a low level of analyte from large samples for accurate analysis. When evaluating analyte extraction or separation, Enviro-Clean® offers nonpolar, polar, ion-exchange and copolymeric phases for application in the environmental laboratory.

Hydrophobic Extraction Columns

Non-polar phases are often referred to as hydrophobic and function by the interactions of the carbon-hydrogen bond of the analyte and the sorbent. C18 is the most widely used of these phases. EPA approved methods for analyzing organics in drinking water specify the C18 hydrophobic phase. This method requires that large sample volumes (liters) be analyzed which utilizes the compound concentration function of the hydrophobic sorbent.

Weight /Vol.	Sorbent	Endcapped	Unendcapped	Sorbent	Endcapped	Unendcapped	Sorbent	Endcapped	Unendcapped	Qty
50 mg/1 mL	C2, Ethyl	EEC021L1	EUC021L1	C4, n-Butyl	EECN41L1	EUCN41L1	C6, Hexyl	EEC061L1	EUC061L1	100
100 mg/1 mL		EEC02111	EUC02111		EECN4111	EUCN4111		EEC06111	EUC06111	100
200 mg/3 mL		EEC02123	EUC02123		EECN4123	EUCN4123		EEC06123	EUC06123	50
500 mg/3 mL		EEC02153	EUC02153		EECN4153	EUCN4153		EEC06153	EUC06153	50
500 mg/6 mL		EEC02156	EUC02156		EECN4156	EUCN4156		EEC06156	EUC06156	50
1000 mg/6 mL		EEC021M6	EUC021M6		EECN41M6	EUCN41M6		EEC061M6	EUC061M6	30
100 mg/10 mL		EEC0211Z	EUC0211Z		EECN411Z	EUCN411Z		EEC0611Z	EUC0611Z	50
200 mg/10 mL		EEC0212Z	EUC0212Z		EECN412Z	EUCN412Z		EEC0612Z	EUC0612Z	50
500 mg/10 mL		EEC0215Z	EUC0215Z		EECN415Z	EUCN415Z		EEC0615Z	EUC0615Z	50
2000 mg/15 mL		EEC0212M15	EUC0212M15		EECN412M15	EUCN412M15		EEC0612M15	EUC0612M15	20
50 mg/1 mL	C3, Propyl	EECN31L1	EUCN31L1	C5, Pentyl	EEC051L1	EUC051L1	C7, Heptyl	EEC071L1	EUC071L1	100
100 mg/1 mL		EECN3111	EUCN3111		EEC0511	EUC05111		EEC07111	EUC07111	100
200 mg/3 mL		EECN3123	EUCN3123		EEC05123	EUC05123		EEC07123	EUC07123	50
500 mg/3 mL		EECN3153	EUCN3153		EEC05153	EUC05153		EEC07153	EUC07153	50
500 mg/6 mL		EECN3156	EUCN3156		EEC05156	EUC05156		EEC07156	EUC07156	50
1000 mg/6 mL		EECN31M6	EUCN31M6		EEC051M6	EUC051M6		EEC071M6	EUC071M6	30
100 mg/10 mL		EECN311Z	EUCN311Z		EEC0511Z	EUC0511Z		EEC0711Z	EUC0711Z	50
200 mg/10 mL		EECN312Z	EUCN312Z		EEC0512Z	EUC0512Z		EEC0712Z	EUC0712Z	50
500 mg/10 mL		EECN315Z	EUCN315Z		EEC0515Z	EUC0515Z		EEC0715Z	EUC0715Z	50
2000 mg/15 mL		EECN312M15	EUCN312M15		EEC0512M15	EUC0512M15		EEC0712M15	EUC0712M15	20

Sample Prep - Solid Phase Extraction

UCT Enviro-Clean® columns



Hydrophobic Extraction Columns

Weight /Vol.	Sorbent	Endcapped	Unendcapped	Sorbent	Endcapped	Unendcapped	Sorbent	Endcapped	Unendcapped	Qty
50 mg/1 mL	C8, Octyl	EEC081L1	EUC081L1	C12, nDodecyl	EEC121L1	EUC121L1	Cyclohexyl	EECYH1L1	EUCYH1L1	100
100 mg/1 mL		EEC08111	EUC08111		EEC12111	EUC12111		EECYH111	EUCYH111	100
200 mg/3 mL		EEC08123	EUC08123		EEC12123	EUC12123		EECYH123	EUCYH123	50
500 mg/3 mL		EEC08153	EUC08153		EEC12153	EUC12153		EECYH153	EUCYH153	50
500 mg/6 mL		EEC08156	EUC08156		EEC12156	EUC12156		EECYH156	EUCYH156	50
1000 mg/6 mL		EEC081M6	EUC081M6		EEC121M6	EUC121M6		EECYH1M6	EUCYH1M6	30
100 mg/10 mL		EEC0811Z	EUC0811Z		EEC1211Z	EUC1211Z		EECYH11Z	EUCYH11Z	50
200 mg/10 mL		EEC0812Z	EUC0812Z		EEC1212Z	EUC1212Z		EECYH12Z	EUCYH12Z	50
500 mg/10 mL	C10, nDecyl	EEC0815Z	EUC0815Z	C18, Octadecyl	EEC1215Z	EUC1215Z	Phenyl	EECYH15Z	EUCYH15Z	50
2000 mg/15 mL		EEC0812M15	EUC0812M15		EEC1212M15	EUC1212M15		EECYH12M15	EUCYH12M15	20
50 mg/1 mL	C10, nDecyl	EEC101L1	EUC101L1	C18, Octadecyl	EEC18111	EUC18111	Phenyl	EEPHY1L1	EUPHY1L1	100
100 mg/1 mL		EEC10111	EUC10111		EEC18123	EUC18123		EEPHY111	EUPHY111	100
200 mg/3 mL		EEC10123	EUC10123		EEC18153	EUC18153		EEPHY123	EUPHY123	50
500 mg/3 mL		EEC10153	EUC10153		EEC18156	EUC18156		EEPHY153	EUPHY153	50
500 mg/6 mL		EEC10156	EUC10156		EEC181M6	EUC181M6		EEPHY156	EUPHY156	50
1000 mg/6 mL		EEC101M6	EUC101M6		EEC1811Z	EUC1811Z		EEPHY1M6	EUPHY1M6	30
100 mg/10 mL		EEC1011Z	EUC1011Z		EEC1812Z	EUC1812Z		EEPHY11Z	EUPHY11Z	50
200 mg/10 mL		EEC1012Z	EUC1012Z		EEC1815Z	EUC1815Z		EEPHY12Z	EUPHY12Z	50
500 mg/10 mL	C10, nDecyl	EEC1015Z	EUC1015Z	C18, Octadecyl	EEC1812M15	EUC1812M15	Phenyl	EEPHY15Z	EUPHY15Z	50
2000 mg/15 mL		EEC1012M15	EUC1012M15					EEPHY12M15	EUPHY12M15	20

Hydrophilic Extraction Columns

Polar or hydrophilic phases function by hydrogen bonding, pi-pi and dipole-dipole interaction. Ion exchange interactions occur between the sorbent and the analyte of opposite charge.

Enviro-Clean® sorbents are available in both cation or anion exchangers exhibiting both weak and strong characteristics.

Copolymeric phases offer a new approach to the environmental analyst by providing very clean extracts and high compound recovery. Dual functionalities, hydrophobic plus ion-exchange or polar allow a higher degree of selectivity than was previously possible. Analytes retained by multiple mechanisms can be washed by disrupting only one mechanism. Careful selection of the solvent strength results in a greater removal of chromatographic contamination.

Weight /Vol.	Qty	P/N
50 mg/1 mL	100	EUCARB1L1
100 mg/1 mL	100	EUCARB111
200 mg/3 mL	50	EUCARB123
500 mg/3 mL	50	EUCARB153
200 mg/6 mL	50	EUCARB126
500 mg/6 mL	50	EUCARB156
1000 mg/6 mL	30	EUCARB1M6
100 mg/10 mL	50	EUCARB11Z
200 mg/10 mL	50	EUCARB12Z
500 mg/10 mL	50	EUCARB15Z
1000 mg/15 mL	20	EUCARB1M15
2000 mg/15 mL	20	EUCARB12M15

Carbon-Graphitized

Application : Carbon supports have been used to isolate extremely polar organic compounds. They work by a hydrophobic mechanism with a high surface area and ion exchange. This interaction can happen in a wide range of polar and non-polar solvents.

Sample Prep - Solid Phase Extraction

UCT Enviro-Clean® columns

Weight /Vol.	Sorbent	P/N	Sorbent	P/N	Sorbent	P/N	Sorbent	P/N	Qty
50 mg/1 mL	Unbonded Silica	EUSIL1L1	Florisil®	EUFLS1L1	Alumina, Acidic	EUALA1L1	Alumina, Basic	EUALB1L1	100
100 mg/1 mL		EUSIL111		EUFLS111		EUALA111		EUALB111	100
200 mg/3 mL		EUSIL123		EUFLS123		EUALA123		EUALB123	50
500 mg/3 mL		EUSIL153		EUFLS153		EUALA153		EUALB153	50
500 mg/6 mL		EUSIL156		EUFLS156		EUALA156		EUALB156	50
1000 mg/6 mL		EUSIL1M6		EUFLS1M6		EUALA1M6		EUALB1M6	30
100 mg/10 mL		EUSIL11Z		EUFLS11Z		EUALA11Z		EUALB11Z	50
200 mg/10 mL		EUSIL12Z		EUFLS12Z		EUALA12Z		EUALB12Z	50
500 mg/10 mL		EUSIL15Z		EUFLS15Z		EUALA15Z		EUALB15Z	50
2000 mg/15 mL		EUSIL12M15		EUFLS12M15		EUALA12M15		EUALB12M15	20
50 mg/1 mL	Alumina, neutral	EUALN1L1	CN, Cyanopropyl Endcapped	EECNP1L1	CN, Cyanopropyl Unendcapped	EUCNP1L1	Diol	EUDOL1L1	100
100 mg/1 mL		EUALN111		EECNP111		EUCNP111		EUDOL111	100
200 mg/3 mL		EUALN123		EECNP123		EUCNP123		EUDOL123	50
500 mg/3 mL		EUALN153		EECNP153		EUCNP153		EUDOL153	50
500 mg/6 mL		EUALN156		EECNP156		EUCNP156		EUDOL156	50
1000 mg/6 mL		EUALN1M6		EECNP1M6		EUCNP1M6		EUDOL1M6	30
100 mg/10 mL		EUALN11Z		EECNP11Z		EUCNP11Z		EUDOL11Z	50
200 mg/10 mL		EUALN12Z		EECNP12Z		EUCNP12Z		EUDOL12Z	50
500 mg/10 mL		EUALN15Z		EECNP15Z		EUCNP15Z		EUDOL15Z	50
2000 mg/15 mL		EUALN12M15		EECNP12M15		EUCNP12M15		EUDOL12M15	20

Ion exchange and mixed mode Extraction Columns

Weight /Vol.	Sorbent	P/N	Sorbent	P/N	Sorbent	P/N	Sorbent	P/N	Qty
50 mg/1 mL	Diethylamino	EUDAX1L1	Aminopropyl	EUNAX1L1	Quaternary Amine with Chloride Counter Ion	EUQAX1L1	Carboxylic Acid	EUCCX1L1	100
100 mg/1 mL		EUDAX111		EUNAX111		EUQAX111		EUCCX111	100
200 mg/3 mL		EUDAX123		EUNAX123		EUQAX123		EUCCX123	50
500 mg/3 mL		EUDAX153		EUNAX153		EUQAX153		EUCCX153	50
500 mg/6 mL		EUDAX156		EUNAX156		EUQAX156		EUCCX156	50
1000 mg/6 mL		EUDAX1M6		EUNAX1M6		EUQAX1M6		EUCCX1M6	30
100 mg/10 mL		EUDAX11Z		EUNAX11Z		EUQAX11Z		EUCCX11Z	50
200 mg/10 mL		EUDAX12Z		EUNAX12Z		EUQAX12Z		EUCCX12Z	50
500 mg/10 mL		EUDAX15Z		EUNAX15Z		EUQAX15Z		EUCCX15Z	50
2000 mg/15 mL		EUDAX12M15		EUNAX12M15		EUQAX12M15		EUCCX12M15	20
50 mg/1 mL	Benzene- -sulfonic Acid	EUBCX1L1	Hydrophobic plus Carboxylic Acid	EUCCX2L1	Hydrophobic plus Benzene- -sulfonic Acid	EUBCX2L1	Hydrophobic plus Quaternary Amine	EUQAX2L1	100
100 mg/1 mL		EUBCX111		EUCCX211		EUBCX211		EUQAX211	100
200 mg/3 mL		EUBCX123		EUCCX223		EUBCX223		EUQAX223	50
500 mg/3 mL		EUBCX153		EUCCX253		EUBCX253		EUQAX253	50
500 mg/6 mL		EUBCX156		EUCCX256		EUBCX256		EUQAX256	50
1000 mg/6 mL		EUBCX1M6		EUCCX2M6		EUBCX2M6		EUQAX2M6	30
100 mg/10 mL		EUBCX11Z		EUCCX21Z		EUBCX21Z		EUQAX21Z	50
200 mg/10 mL		EUBCX12Z		EUCCX22Z		EUBCX22Z		EUQAX22Z	50
500 mg/10 mL		EUBCX15Z		EUCCX25Z		EUBCX25Z		EUQAX25Z	50
2000 mg/15 mL		EUBCX12M15		EUCCX22M15		EUBCX22M15		EUQAX22M15	20