

## DNA Preparation from Paraffin Tissue

Section of Cancer Genomics, Genetics Branch, NCI  
National Institutes of Health

### Reagents

**Chloroform**

Mallinckrodt, Cat. 4440

**EDTA, 0.5 M****Ethanol, absolute****Isoamyl alcohol**

Sigma, Cat. I-3643

**Phenol****Proteinase K****RNase A**

Boehringer, Cat. 109 169

**Sodium acetate, pH 5.2****Sodium chloride, 5 M****Sodium thiocyanate (NaSCN), 1 M**

Sigma, Cat. S 7757, 250 g)

**TE Buffer (Tris-EDTA), pH 7.4****Tween 20****Xylene**

### Preparation

**Chloroform/Isoamyl alcohol 24:1**

Chloroform            24 ml

Isoamyl alcohol       1 ml

**DNA extraction buffer**

1.5 ml 5M NaCL

5.0 ml 0.5M EDTA

0.5 ml Tween 20

Fill up to 100 ml with sterile water.

**Proteinase K (10 mg/ml)**

Dissolve 100 mg Proteinase K in 10 ml TE for 30 min at RT.

Aliquot and store at -20°.

**RNase A (20 mg/ml)**

Dissolve 200 mg RNase A in 5 ml sterile water.

Boil for 15 min.

Cool to room temperature

Aliquot and store at  $-20^{\circ}$ .

**Procedure**

1. Cut 50  $\mu\text{m}$  slices of formalin-fixed and paraffin-embedded tumor samples.  
(Note: Cut 4  $\mu\text{m}$  slices before and after each 50  $\mu\text{m}$  slice for Haematoxylin-Eosin staining to insure that the tissue is still representative.)
2. Incubate in xylene at  $45^{\circ}\text{C}$  for 15 min.
3. Centrifuge 10 min at 14,000 rpm.
4. Pipet off supernatant.
5. Repeat steps 2-4.
6. Add 1 ml 100% ethanol to tissue pellet, vortex, centrifuge for 10 min at 14,000 rpm, and pipet off supernatant.
7. Add 1 ml 90% ethanol to tissue pellet, vortex, centrifuge for 10 min at 14,000 rpm, and pipet off supernatant.
8. Add 1 ml 70% ethanol to tissue pellet, vortex, centrifuge for 10 min at 14,000 rpm, pipet off supernatant, and dry pellet in speed vac.
9. Resuspend pellet in 1 ml NaSCN (1 M) and incubate at  $37^{\circ}\text{C}$  overnight.
10. Centrifuge for 10 min at 14,000 rpm; pipet off supernatant.
11. Resuspend pellet in 400  $\mu\text{l}$  of DNA extraction buffer.
12. Add 5  $\mu\text{l}$  RNase (20 mg/ml) and incubate for 1 hr at  $37^{\circ}\text{C}$  (RNase treatment is optional, paraffin material often does not contain large amounts of RNA).
13. Add 40  $\mu\text{l}$  of Proteinase K (10 mg/ml), vortex briefly, and incubate at  $55^{\circ}\text{C}$  overnight (if tissue is not completely dissolved, add additional proteinase K and continue incubating; tissue should be dissolved).

14. Add 440  $\mu$ l of phenol, shake vigorously by hand for 5 min, and centrifuge for 5 min at 8000 rpm.
15. Pipet supernatant into a new tube, add a solution of 220  $\mu$ l phenol plus 220  $\mu$ l chloroform/isoamyl alcohol (24:1), shake vigorously by hand for 5 min, and centrifuge for 5 min at 8000 rpm.
16. Pipet supernatant into a new tube, add 440  $\mu$ l chloroform/isoamyl alcohol (24:1), shake vigorously by hand for 5 min, and centrifuge for 5 min at 8000 rpm.
17. Pipet supernatant into a new tube (2 ml eppendorf tube), add 1/10 volume of sodium acetate (pH 5.2), add 3 volumes of ice cold 100 % ethanol, and keep tube for 1 hr at  $-80^{\circ}\text{C}$  or overnight at  $-20^{\circ}\text{C}$ .
18. Centrifuge for 30 min at  $4^{\circ}\text{C}$ .
19. Remove and save supernatant (optional, in case you are not sure of the precipitation).
20. Dry pellet in speed vac.
21. Add 20-50  $\mu$ l sterile  $\text{H}_2\text{O}$  (depending on the amount of DNA you expect, which is subjective to experience).
22. Shake gently in thermomixer at  $37^{\circ}\text{C}$  for 2 hours (DNA should be dissolved, but if you have doubts put on a rotating shaker in the cold room overnight).
23. Measure DNA concentration with a spectrophotometer and run around 200 ng on a 1% agarose gel.